

Product data sheet

HT-Fect™ transfection reagent

Catalog Number: TR-101

Components: HT-Fect™ transfection reagent (1 ml)
HT-Fect™ buffer (15 ml)

Storage:

Components: HT-Fect™ transfection reagent (1 ml): -80°C
HT-Fect™ buffer (15 ml): 4°C

Product description

HT-Fect™ transfection reagent is a polymer-based transfection reagent optimized for delivering nucleic acid into mammalian cells. HT-Fect™ has been used in transfecting a broad range of cell types, including hard-to-transfect T cells, fibroblast cells and neuronal cells. HT-Fect™ is recommended to be used with our shRNA plasmids and gRNA plasmids, and routinely used for our high titer lentivirus production.

Transfection Procedure

Note: all quantity and volume are given on per well of 6-well plate, it should be scale up or down for other cell culture dishes.

1. Plate cells

- a. For adherent cells, plate cells one day before the transfection experiment so that the cells will be 60%-80% confluent on the day of transfection.
- b. For suspension cells, suspension culture cells should be in good growth condition before transfection.

2. Prepare the HT-Fect™ transfection reagent

Warm the HT-Fect™ transfection reagent and buffer to room temperature. Mix well before use.

3. Prepare transfection complex

- 1) Solution A: Add 2 ug DNA to a sterile 1.5 ml centrifuge tube and dilute with HT-Fect™ buffer to final volume of 20 µl.
- 2) Solution B: Add 2 ul of HT-Fect™ transfection reagent to a sterile 1.5 ml centrifuge tube and dilute with HT-Fect™ buffer to 20 µl.
- 3) Mix Solution A and B and incubate at room temperature for 15 minutes.

4. Add transfection complex to the cells

- 1) Dilute 40 µl of the above transfection complex by adding 1ml of Opti-MEM or other serum-free growth medium

- 2) a. For adherent cells, Aspirate the growth medium from the wells, and add the above diluted transfection complex to the well.
 - b. For suspension cells, centrifuge the cells at 400 g for 5 minutes, remove the culture medium, and resuspend the cells with the above diluted transfection complex, put the cells into the cell culture plate.
 - 3) Return cells to the incubator.
- 5. Change back to complete cell culture medium 2 hour to overnight post-transfection.**
 - 6. Incubate the cells for 48-96 hours prior to check transgene expression.**